5. Tip-to-Tip Pinch

Place the tip of the thumb on the center of the ergonomic palm bar. Position the finger to be exercised on one of the center buttons. Keeping both tips slightly bent, flex the thumb and the finger being exercised toward the center of the unit at the same time.

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
6. Fingertip Flexion

Place all finger tips on their own button, with the tip of the thumb on the center of the ergonomic base. Keeping all finger tips bent, flex the thumb and all the fingers toward the center of the unit at the same time.

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
PROGRESSIVE EXERCISE SYSTEMS

DIGI-FLEX

exercise each finger individually, or entire hand

Resistance level: per finger / per hand
 Weerstandsniveau: per vinger / geheel
 Niveau de résistance: par doigt / sur l'ensemble
 Widerstandskraft: pro finger / gesamtes gerät
 Livello di resistenza: per singolo dito / intero apparecchio
 Nivel de resistencia: por dedo / sistema completa

10-0751		10-0740		10-0741		10-0742	
tan	xx-light	yellow	x-light	red	light	green	medium
braun	xx-leicht	gelb	x-leicht	rot	leicht	grün	mittel
beige	xx-licht	geel	x-licht	rood	licht	groen	medium
beige	xx-légère	jaune	x-légère	rouge	légère	vert	moyen
marrone	xx-leggero	giallo	x-leggero	rosso	leggero	verde	medio
marfil	xx-ligero	amarillo	x-ligero	rojo	ligero	verde	mediano
	0.75 lbs / 2.5 lbs		1.5 lbs / 5.0 lbs		3.0 lbs / 10.0 lbs		5.0 lbs / 16.0 lbs
	0.34 kgs / 1.1 kgs		0.7 kgs / 2.3 kgs		1.4 kgs / 4.5 kgs		2.3 kgs / 7.3 kgs


10-0743		10-0744		10-0752		10-0753	
blue	heavy	black	x-heavy	silver	xx-heavy	gold	xxx-heavy
blau	stark	schwarz	x-stark	silber	xx-stark	golden	xxx-stark
blauw	zwaar	zwart	x-zwaar	zilver	xx-zwaar	goud	xxx-zwaar
bleu	fort	noir	x-fort	argent	xx-fort	doré	xxx-fort
blu	pesante	nero	x-pesante	argento	xx-pesante	oro	xxx-pesante
azul	pesado	negro	x-pesado	plata	xx-pesado	oro	xxx-pesado
	7.0 lbs / 23.0 lbs		9.0 lbs / 30.0 lbs		11.0 lbs / 36.0 lbs		13.0 lbs / 43.0 lbs
	3.2 kgs / 10.4 kgs		4.1 kgs / 13.6 kgs		5.0 kgs / 16.3 kgs		5.9 kgs / 19.5 kgs

now available:

Cando® Digi-Extend® exerciser for finger extension (a)

Cando® Digi-Squeeze® finger exerciser for hand grasp and finger flexion (b)

Cando® Digi-Extend n'Squeeze® for flexion and extension combination exercises (c)

Manufactured by:

DIGI-FLEX®

a Division of FEI
 P.O. Box 1500
 White Plains, NY 10602 USA

tel: 800-431-2830 • 914-345-9300
 fax: 800-634-5370 • 914-345-9800
 email: info@Fab-Ent.com

Cando®, Digi-Flex® are trademarks of Goldberg
 ©Copyright FEI, 2011 - all rights reserved

Authorized CE Representative:
 RMS UK ltd.
 28 Trinity Road
 Nailsea, Somerset BS48 4NU
 United Kingdom

Digi-Flex® Exercise System:

10-0751, 10-0740, 10-0741, 10-0742, 10-0743, 10-0744, 10-0752, 10-0753

INDICATIONS

- Develops individual finger strength, flexibility and coordination
- Also builds hand and forearm strength

PRECAUTIONS

- **If the individual is being treated for a hand disorder, specific exercises and the intensity and number of repetitions should be approved by a healthcare professional.**
- **Following injury or surgery, consult a healthcare professional before performing exercises.**

- **If unusual swelling, skin discoloration or discomfort occurs, use should be discontinued and a healthcare professional consulted.**

GUIDELINES FOR USE

Exercises should be done at a slow, comfortable pace. Suggested exercises follow.

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1. Gross Grasp

Place the middle of the fingers on the buttons with the finger tips wrapped over the top. Positioning the hook of the ergonomic palm bar over the web space between the thumb and index finger, flex the thumb and all fingers toward the center of the unit as if making a fist.

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
2. Hook Position

Place finger tips on individual buttons and position ergonomic palm bar with the hook over the web space between the thumb and index finger. Keeping finger tips slightly bent, flex all fingers toward the center of the unit.

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
3. Key Pinch

Place the hook of the ergonomic palm bar on the inside of the bent index finger. With the tip of the thumb slightly bent, flex the thumb on the first button.

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
4. Trigger Pinch

Loosely hold the unit vertically with the buttons toward the fingers and the hook of the ergonomic palm bar positioned over the web space between the thumb and index finger. With the middle of the index finger, press the second button as if squeezing a trigger.

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