

# Nebulizer and CPAP Cleaners



## 3B QN1000 Qube Nebulizer Compressor Kit

Sku: QN1000

Power: AC 120V, 60Hz, ETL Classified Listed | Power Consumption: 1.3A | Medication Capacity: 6ml | Particle Size: >75% 0.5~5µm | MMAD: 2.97µm | Sound Level: Less than 60 dB at 1 meter | Average Nebulization Rate: ≥0.2ml/min | Compressor Pressure Range: 30 to 50 PSI | Liter Flow Range: 6.0 ~ 10.0 LPM | Operating/Storage Temperature Range: 10°C to 40°C (50°F to 104°F)



## Philips Respironics 1099966 Innospire Essence w/ SideStream Nebulizer

Sku: 1099966

Comes with SideStream Disposable Nebulizer | Maximum flow rate: 9.4 LPM | Maximum pressure: 44 psi | Mode of operation: 30 min. on/30 min. off min. | Nebulization capacity: 8 ml | Operating flow rate: 7 LPM @ 10 psi LMP/psi | Size: 4.2 x 6.5 x 6.5 | Treatment time: 5 to 7 minutes | Weight: 3.3 lbs



## Philips Respironics 1100312 Innospire Essence w/ SideStream Nebulizers

Sku: 1100312

Comes with SideStream Disposable & Reusable Nebulizer | Electrical rating/power: 120 V / 60 Hz, 1.7 Amp V/Hz/Amp | Maximum flow rate: 9.4 LPM | Maximum pressure: 44 psi | Mode of operation: 30 min. on/30 min. off min. | Nebulization capacity: 8 ml | Operating flow rate: 7 LPM @ 10 psi LMP/psi | Treatment time: 5 to 7 minutes | Weight: 3.3 lbs



## Flyp PE1200M Portable Vibrating Mesh Nebulizer

Sku: PE1200M

Average Treatment Time: 7 minutes | Dimensions: 1.1" x 4.7" x 2.1" | Disk Characteristics: Stainless steel 1000 precision laser cut 4 micron holes | Method of Nebulization: HypersonIQ Vibrating Mesh Disk | Power Supply: Lithium Ion micro-USB rechargeable battery | Treatments Per Charge: 8/10/2019 | Weight: 3.6 oz



## 3B Medical Lumin CPAP Mask and Accessories Cleaner

Sku: LM3000

Safely Disinfects Equipment | Effective UV Light System | Quick, Five-Minute Cycle | Non-CPAP Item Sanitization | Easy-to-Use Design | No Scheduled Replacement Items



## Responsive Respiratory Purify O3 CPAP/BiPap Cleaner Sanitizer

Sku: 190-6000

Includes: Purify O3 Sanitizer, dual purpose sanitizing/travel case, tubing adapter, USB cable, wall adapter, bilingual instruction manual | FDA Registered | FSA/HSA Eligible | May be used for disinfection purposes with other small reusable medical equipment including nebulizers, breast pumps, etc. | No extra supplies required - kit includes all you need | Powerful, hospital grade sanitation cycle utilizes ozone to kill 99.99% of germs and bacteria

\*Inventory levels are fluid. Please contact us for current availability

TO PLACE AN ORDER  
 Health Mobius  
 630-325-5150. | [cs2@healthmobius.com](mailto:cs2@healthmobius.com)  
[wholesale.healthmobius.com](http://wholesale.healthmobius.com)

Available on the MHA, Navigator and AHP Contracts

